

# Restwell Medical Supply - CPAP Mask Comfort Tips

Restwell Medical Supply - Patient Comfort Sheet

For: ResMed AirSense 11 users with the AirFit P30i nasal pillow

Getting a better seal on the road

- Fit the pillows sitting up, then lie down and adjust. Gravity changes the seal.
- A worn cushion leaks. Rotate on schedule, keep a spare in the sleeper.
- Wipe the cushion each morning. Diesel dust in a truck stop lot is real.

Sleeper cab notes

- Keep the hose off the floor. A hose clip on the bunk rail helps.
- Cold air stiffens the cushion. Let it warm a minute before you seal up.

Compliance is what keeps the DOT card. Your data uploads through the myAir app.

Questions on supplies go through your Restwell account. Talk to Dr. Cho about therapy.